

# ENTHIUS

impression  
materials



## ENTHIUS

is an **elastomeric vinyl polysiloxane impression material** that provides the clinician with an **accurate and reliable impression** with exceptional tear strength and stability as well as exceptional accuracy and long term dimensional stability.

- Enthus impression materials feature **high fluidity** and **hydrophilicity** with high contrast colors.
- Enthus is a thixotropic material that exhibits minimal increase in viscosity, remaining **flowable during working time**. As soon as the tray is placed in the mouth Enthus' viscosity will increase immediately.
- Enthus offers a unique combination of **high tear strength** and very **low contact angle**.
- Gluten-Free
- Made in USA
- Available in Regular & Fast

# Impression materials **ENTHUS**

AVAILABLE	COLOR	WORKING TIME (minutes)		INTRA-ORAL SET TIME (minutes)		SHORE A HARDNESS
		REGULAR SET	FAST SET	REGULAR SET	FAST SET	
Light Body	Yellow	2:30	1:00	1:00	1:00	40
Medium Body	Buff	2:30	1:00	1:00	1:00	50
Heavy Body	Rose	2:30	1:00	1:00	1:00	60
Monophase	Purple	2:30	1:00	1:00	1:00	60
Bite Registration	Light Blue	1:00	0:30	1:00	1:00	80
Putty	Dark Blue	2:00	1:00	2:00	2:00	60

## presentations

### Putty: 2 x 300ml jars & 2 scoops

- Item #38-00001  
VPS Putty, Regular Set  
Spearmint
- Item #38-00002  
VPS Putty, Fast Set  
Spearmint



### Bite Registration: 2 x 50ml cartridges, 6 mixing tips

- Item #38-00011  
Bite Registration, Regular Set  
Spearmint
- Item #38-00012  
Bite Registration, Fast Set  
Spearmint



### VPS: 4 x 50ml cartridges, Regular Set

- Item #38-00003  
VPS Light Body, Spearmint
- Item #38-00004  
VPS Medium Body, Spearmint
- Item #38-00005  
VPS Heavy Body, Spearmint
- Item #38-00006  
VPS Monophase, Spearmint



### VPS: Accessories:

- Item #38-00013 - Dispensing Gun 1:1
- Item #38-00014 - Mixing Tips Purple - 7.5mm, Bag x 48
- Item #38-00015 - Mixing Tips Teal - 6.5mm, Bag x 48
- Item #38-00016 - Mixing Tips Pink - 5.4mm, Bag x 48
- Item #38-00017 - Mixing Tips Yellow - 4.2mm, Bag x 48
- Item #38-00018 - Intra-oral Tips Clear, Bag x 100
- Item #38-00019 - Intra-oral Tips Yellow, Bag x 100

### VPS: 4 x 50ml cartridges, Fast Set

- Item #38-00007  
VPS Light Body, Spearmint
- Item #38-00008  
VPS Medium Body, Spearmint
- Item #38-00009  
VPS Heavy Body, Spearmint
- Item #38-00010  
VPS Monophase, Spearmint



Distributed by: